

CHART 2020

Your practice logo goes here



THE PREMIER VISION ACUITY AND OCULAR PERFORMANCE TESTING SOLUTION

Eye Care USA
Tel: +1 (555) 555 5555

John Doe
patient name

14
age

2017/08/03 | 14:50
date | time

Patient Information

9.0.20
version

ETDRS
method

2000/01/01
last calibration date

4.8 Meters
distance


0 sec
inter target pause

Sloan letters
optotype


Test Conditions

	VA 100% Correct	VA 0% Correct	# Letters Correct (Nc)	# Letters Incorrect (Ni)	T = Ni X 0.02	Base LogMar	LogMAR VA = Base LogMAR + T
OD ●	100	32	28	10	0.20	0.3	0.5
OS ■	125	32	23	17	0.34	0.3	0.64
OU ▲	50	20	39	11	0.22	0.1	0.32

Summary results



Graphical Overview - Patient VA Ranges



letters correct for each VA

Notes:

Test was performed while patient was dilated

Add customized notes and submit to EHR